Антитела

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волгоград (844)278-03-48 Вологда (8172)26-4159 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89 Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48

Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81 Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73

Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47 Ростов-на-Дону (863)308-18-15 Рязань (4912)46-61-64 Самара (846)206-03-16 Саранск (8342)22-96-24 Санкт-Петербург (812)309-46-40 Саратов (845)249-38-78 Севастополь (8692)22-31-93 Симферополь (3652)67-13-56 Смоленск (4812)29-41-54 Сочи (862)225-72-31 Ставрополь (8652)20-65-13 Сургут (3462)77-98-35

Сыктывкар (8212)25-95-17 Тамбов (4752)50-40-97 Тверь (4822)63-31-35 Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удэ (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

Адрес: https://tci.nt-rt.ru/ || эл.почта: tic@nt-rt.ru

Antibodies and Related Reagents

Antibodies

Anti-Glycolipid Antibodies, Anti-Sulfated Glycan Antibodies	3
Anti-Glycosaminoglycan Antibodies, Anti-Blood Group Antigen Antibodies, Antigen Sugar-conjugated Proteins, Synthetic Carbohydrate Chains	4
Anti-NeuGc Polyclonal Antibodies	5
Anti-αGal Polyclonal Antibodies	6
Anti-Protein A Antibodies, Protein A Products	7
Anti-Tag Antibodies, Anti-Cell Marker Antibodies	8
Anti-Endo-M Antibodies, Anti-Influenza Virus Antibodies	9
Secondary Antibodies and Other Antibodies, Streptavidins	10
Fluorescent Labeled Secondary Antibodies and Fluorescent Cell Stains	11

Products for Labeling

Europium Fluorophore DTBTA-Eu ³⁺ -labeled Proteins	12
Fluorescent Organosilica Particles	13
Peroxidase (HRP) Labeling Reagents	14
Protein-maleimide Conjugates for Thiol-maleimide Crosslinking	16

Products for Cell Biology

Cell Proliferation Assay Reagents, Cell Staining Reagents	17
Extraction Buffers for Cells	18
Enzyme Substrates	
Peroxidase Substrates	19
Alkaline Phosphatase Substrates	22

Other Products

Protein Staining Reagent	23
Protein Determination Reagents	24

Carbohydrate chains are called the third life chain following the protein and the nucleic acid and are one of the most important issues in the post genome research. Most carbohydrate chains attach to lipids or proteins and occur in the form of glycoproteins or glycolipids (*N*-glycan, *O*-glycan, proteoglycans and others). Carbohydrate chains are known to be expressed on brain, nerve, cancer, and endothelial cells. Some carbohydrate chains are known to relate to diseases (e.g., cancer, Alzheimer's disease, Guillain-Barré syndrome, Lysosome syndrome such as Fabry disease, gangliosidosis), differentiation and development (iPS/ES cells). Seasonal influenza viruses, annual epidemics that peak during winter, cause infection via cell-surface glycans. Anti-influenza virus drugs are structural mimics of sialic acid, because neuraminidase is a sialic acid hydrolase that is essential for the release of progeny virus particles from the surface of an infected cell.

Anti-carbohydrate antibodies can recognize glycolipids or glycoproteins. These antibodies can be used for immunohistochemistry, cell-staining, inhibition assay for cell adhesion, flow cytometry, ELISA, TLC-immunostaining and other methods.

Anti-Glycolipid Antibodies

Product Name	lsotype	Size	Product Number
Anti-GM1 Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2505]
Anti-GM ₂ Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2576]
Anti-GM3 Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2582]
Anti-GD _{1a} Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2507]
Anti-GD _{1b} Monoclonal Antibody	Mouse IgG3	0.1mg/vial	[A2508]
Anti-GD ₂ Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A3338]
Anti-GD ₃ Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2580]
Anti-GT _{1a} Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2702]
Anti-GT _{1b} Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2732]
Anti-GQ _{1b} Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2662]
Anti-GalNAc-GD _{1a} Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2701]
Anti-Gb₃ Monoclonal Antibody	Mouse IgG2b	0.1mg/vial	[A2506]
Anti-Gb3 Monoclonal Antibody Biotin Conjugate	Mouse IgG2b	0.1mg/vial	[A2822]
Anti-SGPG (HNK-1) Monoclonal Antibody	Mouse IgG2a	0.1mg/vial	[A2706]

Anti-Sulfated Glycan Antibodies

Product Name	lsotype	Size	Product Number
Anti-6-sulfo LacNAc Monoclonal Antibody (AG105)	Mouse IgM	0.1mg/vial	[A3251]
Anti-6,6'-disulfo LacNAc Monoclonal Antibody (L4L4-8)	Mouse IgM	0.1mg/vial	[A3252]
Anti-Sialyl 6,6'-disulfo LacNAc Monoclonal Antibody (G270-16)	Mouse IgM	0.1mg/vial	[A3253]
Anti- Sialyl 6-sulfo Lewis X Monoclonal Antibody (G152)	Mouse IgM	0.1mg/vial	[A3399]

Anti-Glycosaminoglycan Antibodies			
Product Name	lsotype	Size	Product Number
Anti-Chondroitin Sulfate A Monoclonal Antibody (LY111)	Mouse IgM	0.1mg/vial	[A3143]
Anti-Chondroitin Sulfate D Monoclonal Antibody (MO-225)	Mouse IgM	0.1mg/vial	[A2872]
Anti-Keratan Sulfate Monoclonal Antibody (R-10G)	Mouse IgG1	0.1mg/vial	[A2968]
Anti-Perlecan Monoclonal Antibody (HK-102)	Rat lgG2a	0.1mg/vial	[A3342]

Anti-Blood Group Antigen Antibodies

Product Name	lsotype	Size	Product Number
Anti-Lewis X Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2578]
Anti-Lewis Y Monoclonal Antibody	Mouse IgG3	0.1mg/vial	[A2510]
Anti-Sialyl Lewis A Monoclonal Antibody (1H4)	Mouse IgG3	0.1mg/vial	[A2584]
Anti-Sialyl Lewis A Monoclonal Antibody (2D3)	Mouse IgM	0.1mg/vial	[A2509]
Anti-Sialyl Lewis X Monoclonal Antibody	Mouse IgM	0.1mg/vial	[A2849]

Antigen Sugar-conjugated Proteins

TCI offers carbohydrate-conjugated human serum albumin (HSA) which is manufactured using high-purity synthesized carbohydrates. Several sugar-conjugates are available, and it is also possible to manufacture the sugar-conjugates according to customer specifications. For more details on the products and contracts, please contact us.

HSA-GD ₃
HSA-Lewis
HSA-Sialyl Lewis X
HSA-GM ₁ Pentasaccharide
HSA-Globo-H
HSA-L1-L1

0.1mg/vial [H1718] 0.1mg/vial [H1719] 0.1mg/vial [H1730] 0.1mg/vial [H1767] 0.1mg/vial [H1794]

0.1mg/vial [H1782]



Feel free to contact us. Using advanced proprietary technologies, we synthesize a wide range of sugar chains for daily research.

Anti-NeuGc Polyclonal Antibodies

N-Acetylneuraminic Acid (NeuAc) and *N*-Glycolylneuraminic Acid (NeuGc) are the two major forms of sialic acid found in mammals. Humans are unable to synthesize Neu5Gc due to a mutation in the gene encoding the enzyme responsible for Neu5Gc synthesis. Humans naturally possess antibodies against Neu5Gc glycan structures, and this is responsible for the immunogenicity of therapeutic proteins containing Neu5Gc glycan epitopes. Therefore, a method for the detection of Neu5Gc is required.

Anti-NeuGc Polyclonal Antibody Anti-NeuGc Polyclonal Antibody Biotin Conjugate Anti-NeuGc Polyclonal Antibody FITC Conjugate Anti-NeuGc Polyclonal Antibody R-PE Conjugate Anti-NeuGc Polyclonal Antibody HRP Conjugate 0.05mg/vial [A3240] 0.05mg/vial [A3294] 0.05mg/vial [A3295] 0.05mg/vial [A3360] 0.05mg/vial [A3397]



Detection of NeuGc in miniature pig granulocytes by flow cytometry



Granulocytes were collected by hemolyzing the blood of miniature pigs. The granulocytes were incubated (4 °C, 20 min) with isotype control (black line) or anti-NeuGc polyclonal antibody R-PE conjugate [A3360] (red line) adjusted to 0.1 μ g/mL. Afterward, it was measured using a flow cytometer.

Anti-αGal Polyclonal Antibodies

Anti- α Gal antibody exists as a natural antibody in humans. Binding of this antibody to α Gal antigens (α Gal epitope) expressed on porcine xenograft surfaces are a major factor for determining engraft survival. Recently, it has been observed that therapeutic antibodies and cell processing material for reproductive medicine contain the α Gal epitope, which indicates the importance of rapid detection of α Gal epitope.

Anti-¤Gal Polyclonal Antibody (Chicken)	0.05mg/vial [A3123]
Anti-αGal Polyclonal Antibody Biotin Conjugate	0.05mg/vial [A3144]
Anti-αGal Chicken Polyclonal Antibody HRP Conjugate	0.05mg/vial [A3195]
Anti-¤Gal Polyclonal Antibody FITC Conjugate	0.05mg/vial [A3337]
Anti-αGal Polyclonal Antibody R-PE Conjugate	0.05mg/vial [A3354]



Anti-αGal antibody shows the same high specificity compared with an anti-αGal monoclonal antibody



Glycoconjugates coated on ELISA plates. Results following epitope and anti-αGal antibodies incubation. Primary antibodies were detected using appropriate secondary antibodies.

Detection of αGal in miniature pig granulocytes by flow cytometry



Granulocytes were collected by hemolyzing the blood of miniature pigs. The granulocytes were incubated (4 °C, 20 min) with isotype control (black line) or anti- α Gal polyclonal antibody R-PE conjugate [A3354] (red line) adjusted to 0.1 µg/mL. Afterward, it was measured using a flow cytometer.

Anti-Protein A Antibodies Anti-Protein A Chicken Polyclonal Antibody 0.1mg/vial [A3044] Anti-Protein A Chicken Polyclonal Antibody Biotin Conjugate 0.05mg/vial [A3045] Anti-Protein A Chicken Polyclonal Antibody HRP Conjugate 0.05mg/vial [A3187] High-sensitive detection of Protein A by sandwich-ELISA 1) Dilute anti-Protein A antibody [A3044] with sodium carbonate buffer (pH 8.5), and coat on an ELISA plate. 2) Block with 1% BSA / TBS-T for 2 hours. 3) After washing 3 times with TBS-T, add the sample to each well and incubate Example of calibration curve for 30 minutes. 1.2 4) After washing 3 times with TBS-T, add 1 µg/mL of anti-Protein A antibody 1.0 $R^2 = 0.9997$ biotin conjugate [A3045] to each well and incubate for 30 minutes. Е 0.8 450 5) After washing 3 times with TBS-T, add SA-HRP [S0972] to each well and 0.6 incubate for 30 minutes. Abs. 0.4 6) After washing 3 times with TBS-T, add TMB solution and incubate for 30 minutes. 0.2 7) Stop the reaction by adding 1 N HCI, and measure the absorbance at 450nm. 0.0 0.5 1.5 [ng/mL]

Protein A Products

5mg/vial [P2366]
1mg/vial [P2407]
0.2mg/vial [P2466]
2mL/vial [P2461]



Protein A is a type I membrane protein produced by several strains of *Staphylococcus aureus*. It has high-affinity binding sites for IgGs obtained from various species such as humans, rabbit, mouse, and bovine. Protein A supported by agarose resin is prepared using a covalent coupling method and can be applied to the purification of IgGs. By using **P2461**, human IgGs can be eluted under milder conditions (such as at pH 4.0) compared to using other resins with conventional eluting protocols.

Anti-Tag Antibodies

Anti-DYKDDDDK Antibody Mouse Anti-DYKDDDDK Monoclonal Antibody 0.1mg/vial [M3389] Mouse Anti-DYKDDDDK Monoclonal Antibody Biotin Conjugate 0.05mg/vial [M3400] Anti-HHHHHH (6xHis) Antibody Anti-6xHis Monoclonal Antibody (6A12) 0.1mg/vial [A2957] Immunogen : HHHHHH (6xHis) Isotype : Mouse IgG1 Anti-6xHis Monoclonal Antibody (6A12) Biotin Conjugate 0.05mg/vial [A3010] Anti-6xHis Monoclonal Antibody (6A12) HRP Conjugate 0.05mg/vial [A3075] Anti-Glutathione S-Transferase (GST) Antibody Anti-GST Monoclonal Antibody 0.1mg/vial [A3175] Immunogen : Glutathione S-transferase (GST) Isotype : Mouse IgG2a **Anti-GST Monoclonal Antibody Biotin Conjugate** 0.05mg/vial [A3226] **GST Detection by A3175** DYKDDDDK Tag Detection by M3400 Western blotting Western blotting (human ITGAX fused with DYKDDDDK tag : 0.1mg/lane) (GST: 0.1 µg/lane) kDa kDa 45 Lane M : Molecular weight marker 229.3 Lane 1 : human ITGAX fused with 32 DYKDDDDK tag (reduction) 136.4 2 Μ 1 94.6 Lane M : Molecular weight marker

Anti-Cell Marker Antibodies

71.3

Μ

Mouse Anti-NeuN Monoclonal Antibody

Lane 1 : GST (nonreduction) Lane 2 : GST (reduction)

0.1mg/vial [M3586]

NeuN (RNA binding protein fox-1 homolog 3) is a nuclear protein mainly expressed in postmitotic neurons. Anti-NeuN antibodies are useful markers of mature neurons and widely used in embryology and neuroscience.

Immunofluorescence of Mouse Tissue Section Stained Using M3586Primary antibody:Mouse Anti-NeuN Monoclonal Antibody [M3586]Secondary antibody:Goat Anti-Mouse IgG, Fab Fragment
Cyanine 3 Conjugate [G0598]4 μg of M3586 and 3 μg of G0598 were mixed and incubated for
1.5 hours at 37°C. The mixture was diluted 500 times, added to a
mouse brain section, and incubated overnight at room temperature
with shaking. After washing, sections were observed via a
fluorescence microscope.

Anti-Endo-M Antibodies	
Anti-Endo-M Polyclonal Antibody Immunogen : <i>endo</i> -β-N-Acetylglucosaminidase (Endo-M) Isotype : Rabbit IgG	0.2mg/vial [A2958]
Anti-Endo-M Polyclonal Antibody Biotin Conjugate	0.1mg/vial [A2959]
Enzymes which Transfers the Intact Oligosaccharides	
endo-β-N-Acetylglucosaminidase (=Endo-M) Recombinant: from <i>Mucor hiemalis</i> expressed in <i>Candida boidinii</i>	100m units/vial [A1651]
Glycosynthase (Endo-M-N175Q) Recombinant: from <i>Mucor hiemalis</i> expressed in <i>Escherichia coli</i>	100m units/vial [G0365]
Endo-M-W251N Recombinant: from <i>Mucor hiemalis</i> expressed in <i>Escherichia coli</i>	500m units/vial [E1339]

Anti-Influenza Virus Antibodies

Anti-Influenza A Virus Neuraminidase N1 Monoclonal Antibody			
Immunogen : Influenza A/Brijing/262/95	Clone name : 2-3B	lsotype : Mouse IgĞ1	0.2mL [A2407]
Anti-Influenza A Virus Hemaggl	utinin H3 Mond	clonal Antibody	
Immunogen : Influenza A/Sydney/5/97	Clone name : 1G8	lsotype : Mouse IgĠ3	0.2mL [10779]
Anti-Influenza A Virus Neurami	nidase N2 Mono	oclonal Antibody	
Immunogen : Influenza A/Sydney/5/97	Clone name : 1-4B	Isotype : Mouse IgG1	0.2mL [A2380]
Anti-Influenza A Virus Nucleoprotein Monoclonal Antibody			
Immunogen : Influenza A/Beijing/262/95	Clone name : 17	Isotype : Mouse IgG2a	0.2mL [A2406]

Secondary Antibodies and Other Antibodies

Anti-Mouse IgG Goat Anti-Mouse IgG 1mg/vial [G0386] **Goat Anti-Mouse IgG Biotin Conjugate** 0.1mg/vial [G0387] Goat Anti-Mouse IgG HRP Conjugate 0.1mg/vial [G0407] **Goat Anti-Mouse IgG FITC Conjugate** 0.1mg/vial [G0406] Goat Anti-Mouse IgG R-PE Conjugate 0.1mg/vial [G0569] Goat Anti-Mouse IgG, Fab Fragment Cyanine 3 Conjugate 0.05mg/vial [G0598] Anti-Mouse IgM Goat Anti-Mouse IgM 1mg/vial [G0408] Goat Anti-Mouse IgM Biotin Conjugate 0.1mg/vial [G0432] Goat Anti-Mouse IgM HRP Conjugate 0.1mg/vial [G0417] Goat Anti-Mouse IgM FITC Conjugate 0.1mg/vial [G0453] Anti-Rabbit IgG Goat Anti-Rabbit IgG 1mg/vial [G0388] Goat Anti-Rabbit IgG Biotin Conjugate * 0.1mg/vial [G0597] Goat Anti-Rabbit IgG HRP Conjugate 0.1mg/vial [G0418] Goat Anti-Rabbit IgG FITC Conjugate 0.1mg/vial [G0452] Goat Anti-Rabbit IgG R-PE Conjugate 0.1mg/vial [G0577] Anti-Chicken IgY **Sheep Anti-Chicken IqY** 1mg/vial [S0998] Sheep Anti-Chicken IgY Biotin Conjugate 0.1mg/vial [H1619] Sheep Anti-Chicken IgY HRP Conjugate 0.1mg/vial [S0999] Anti-HRP Antibody **Anti-HRP Rabbit Polyclonal Antibody** 0.2mL [A2250] Immunogen : Horseradish Peroxidase Isotype: Rabbit IgG Anti-Human lgG Anti-Human IgG Fc C-terminus Monoclonal Antibody 0.1mg/vial [A3277] Immunogen : Synthetic peptide corresponding to human IgG Fc C terminus Isotype: MouseIgG1 Mouse Anti-Human IgG Fc 0.1mg/vial [M2977] Mouse Anti-Human IgG Fc Biotin Conjugate 0.1mg/vial [M3053] *G0597 is the successor to Anti-Rabbit IgG Biotin Conjugate (Product Number: G0389). Please use G0597 alternatively if you have used G0389. **Streptavidins** Streptavidin from Streptomyces avidinii 1mg/vial [S0951]

Streptavidin from Streptomyces avidin Streptavidin HRP Conjugate Streptavidin FITC Conjugate Streptavidin DTBTA-Eu³⁺ Conjugate Streptavidin R-PE Conjugate Streptavidin Maleimide Conjugate 1mg/vial [S0951] 0.1mg/vial [S0972] 0.1mg/vial [S0966] 0.1mg/vial [S0993] 0.1mg/vial [T3885] 0.5mg/vial [T3531]

Fluorescent Labeled Secondary Antibodies and Fluorescent Cell Stains

Applications



- (A) The HeLa cells were incubated with properly diluted primary antibody (Mouse Anti α-Tubulin IgG) and were further incubated with Goat Anti-Mouse IgG Biotin Conjugate [60387] and Streptavidin FITC Conjugate [S0966] (green fluorescence). And then the nuclei was stained with DAPI 2HCI [A2412] (blue fluorescence). (Laser Scanning Microscope: Olympus FLUOVIEW FV 3000)
- (B) The nuclei of HeLa cells was stained with Bisbenzimide H 33258 [H1343] (blue fluorescence). α-Tubulin was stained with anti-α-tubulin antibody and Goat Anti-Mouse IgG Biotin Conjugate [G0387] and Streptavidin R-PE Conjugate [T3885] (red fluorescence). Mitochondria was stained with primary antibody and Goat Anti-Rabbit IgG FITC Conjugate [G0452] (green fluorescence)**. (Laser Scanning Microscope: Olympus FLUOVIEW FV 3000)



(C) The HeLa cells were incubated with Mouse Anti-CD9 Antibody (red line) or Mouse IgG2ak isotype control (black line). Subsequently, both were stained with Goat Anti-Mouse IgG Biotin Conjugate [G0387] and Streptavidin R-PE Conjugate [T3885]. (Flow cytometer: Sysmex RF-500)

**Please refer to our product page for staining procedure. R-PE/FITC-labeled anti-Mouse IgG or anti-Rabbit IgG antibodies and streptavidins can be used for fluorescence immunostaining and flow cytometry.

Goat Anti-Mouse IgG FITC Conjugate	(Green Fluorescence	0.1mg/vial	[G0406]
Goat Anti-Mouse IgM FITC Conjugate	(Green Fluorescence	0.1mg/vial	[G0453]
Goat Anti-Rabbit IgG FITC Conjugate	(Green Fluorescence	0.1mg/vial	[G0452]
Streptavidin FITC Conjugate	(Green Fluorescence	0.1mg/vial	[S0966]
Goat Anti-Mouse IgG R-PE Conjugate	(Red Fluorescence)	0.1mg/vial	[G0569]
Goat Anti-Mouse IgG, Fab Fragment Cyanine 3 Conjugate	(Red Fluorescence)	0.05mg/vial	[G0598]
Goat Anti-Rabbit IgG R-PE Conjugate	(Red Fluorescence)	0.1mg/vial	[G0577]
Streptavidin R-PE Conjugate	(Red Fluorescence)	0.1mg/vial	[T3885]
Goat Anti-Mouse IgG DTBTA-Eu ³⁺ Conjugate	(Red Fluorescence)	0.1mg/vial	[G0505]
Goat Anti-Rabbit IgG DTBTA-Eu ³⁺ Conjugate	(Red Fluorescence)	0.1mg/vial	[G0506]
Streptavidin DTBTA-Eu ³⁺ Conjugate	(Red Fluorescence)	0.1mg/vial	[S0993]
DAPI 2HCI	(Blue Fluorescence)	5mg	[A2412]
DAPI 2HCI (1mg/mL in Water)	(Blue Fluorescence)	0.2mL x 5vial	[D5888]
Bisbenzimide H 33258 Hydrate	(Blue Fluorescence)	25mg	[H1343]
Bisbenzimide H 33258 (1mg/mL in Water)	(Blue Fluorescence)	0.2mL x 5vial	[B6236]

*Some products are unavilable in the Americas and China. *The high-sensitivity detection of DTBTA-Eu³⁺ labeled probes requires time-resolved fluorometry.

Europium Fluorophore DTBTA-Eu³⁺-labeled Proteins

Highly-sensitive Detection Probes for Time-resolved Fluorometry

Goat Anti-Mouse IgG DTBTA-Eu³⁺ Conjugate Goat Anti-Rabbit IgG DTBTA-Eu³⁺ Conjugate Streptavidin DTBTA-Eu³⁺ Conjugate 0.1mg/vial [G0505] 0.1mg/vial [G0506] 0.1mg/vial [S0993]

Advantages

No cross talk of excitation light

+ Excitation wavelength Ex_{max} : 335 nm

- Emission wavelength Emmax : 616 nm
- Sharpened emission spectrum

Large Stokes shift (the difference in wavelength between positions of the band maxima of the absorption and emission spectra)

Stable fluorescence in various aqueous buffers

Available in Tris, TE, PBS, etc., for wide use

Long fluorescent life time (τ = 1.02 ms)

Time-resolved fluorometric measurement can remove background fluorescence from the sample matrix and often gives detectability better than one order of magnitude compared to those of conventional fluorometric assays.

Comparison of secondary antibody conjugated to DTBTA-Eu³⁺ or FITC

Time-resolved fluorometric measurement can remove background fluorescence! To obtain a high SN ratio



<Assay condition>

Dilute the Mouse IgG to each concentration. Coat 96-well plates with diluted Mouse IgG. Block the plates with BSA/TBST. Incubate with Goat Anti-Mouse IgG Conjugates labeled by DTBTA-Eu³⁺ or FITC at 2.5 µg/mL. After incubation, measure the fluorescence intensity on a plate reader. DTBTA-Eu³⁺; excitation=340 nm, emission=620 nm. Lag Time : 450 µsec FITC; excitation=485 nm, emission=520 nm.

Anti-DTBTA-Eu³⁺ Antibody Anti-DTBTA-Eu³⁺ Rabbit Polyclonal Antibody Anti-DTBTA-Eu³⁺ Rabbit Antiserum

DTBTA-Eu³⁺ Labeling Reagent ATBTA-Eu³⁺

0.5mL [A2239] 0.5mL [A2181]

10mg [A2083]

ositions a)



Fluorescent Organosilica Particles

Organosilica FITC (100nm Diam.) Organosilica Rhodamine B (100nm Diam.)

2mg [00561] 2mg [00573]





Fluorescent image of 00561

SEM image

Advantages

- Wavelength : Ex_{max} 492 nm, Em_{max} 523 nm (O0561) Ex_{max} 556 nm, Em_{max} 579 nm (O0573)
- Surface Functionalization : Thiol group (-SH)
- Superior in fluorescence intensity to the conventional FITC or rhodamine B.
- The diameter of these products are 100 nm and these products are suitable for the detection of biomolecules.



Peroxidase (HRP) Labeling Reagents

Horseradish Peroxidase Maleimide Conjugate (0.5mg×3) Horseradish Peroxidase NHS Ester Conjugate (0.2mg×3) 1set [H1621] 1set [H1746]

Advantages

- H1746 contains an *N*-hydroxysuccinimidyl ester (NHS) moiety and can be used to readily label proteins and peptides that have an amino group (-NH₂).
- H1621 can be used for the conjugation to free thiol-containing proteins and peptides due to its thiol-reactive maleimide group.
- Each protein conjugate is packaged for single use purposes and thus does not require weighing prior to use.



Application : HRP-labelling of an antibody with H1621

In case of antibodies without free thiol (SH, sulfhydryl) groups, disulfide moieties in proteins can be reduced by a reductant such as DTT [D3647] or 2-MEA [A0296] to reveal free thiols. Furthermore, thiol group can be introduced to primary amines by adding SATA [S0431], SATP [S0859] or 2-Iminothiolane.



Example protocol for antibody conjugation starts from a reduction of native disulfide bonds in the Goat Anti-Mouse IgG, followed by labeling with the HRP using H1621. For more information, see the product detail page of H1621 on TCI website.

Protocol

- 1) Add DTT to a final concentration equal to 3 mole equivalents per mole equivalent of antibody present.
- 2) Incubate for 90 minutes at 37 °C.
- 3) Purify the reduced IgG by gel filtration or ultrafiltration, dialysis.
- 4) Add equal amount of H1621 (by weight) to a purified antibody and incubate for 2 hours at room temperature (25 °C).



Goat Anti-Mouse IgG labeled with the HRP using H1621 was tested by ELISA for detection of a Mouse IgG coated on a plate. Mouse IgG could be detected sufficiently even if the labeled antibody was diluted to 5 ng/mL or more.

Related Products

(Reducing agents for protein disulfide)

DTT (= DL-Dithiothreitol) 2-MEA (= 2-Aminoethanethiol Hydrochloride) 2-Mercaptoethanol Tris(2-carboxyethyl)phosphine Hydrochloride

Reagents for introduction of thiol group

SATA (= N-Succinimidyl S-Acetylthioglycolate) SATP (= N-Succinimidyl 3-(Acetylthio)propionate) 1g / 5g [D3647] 25g / 100g / 500g [A0296] 5g / 25g [M1948] 1g / 5g / 25g [T1656]

> 1g / 5g [**S0431**] 100mg [**S0859**]

Protein-maleimide Conjugates for Thiol-maleimide Crosslinking

Bovine Serum Albumin Maleimide Conjugate (1mg×3) Horseradish Peroxidase Maleimide Conjugate (0.5mg×3) Streptavidin Maleimide Conjugate (0.5mg×1) 1set [**B5944**] 1set [**H1621**] 1vial [**T3531**]

Advantages

- Each product containing a thiol-reactive maleimide group can be used for the conjugation to proteins and peptides containing free thiols.
- Each protein conjugate is packaged for single use purposes and thus does not require weighing prior to use.





KIPA Butter (Read	v-to-use) [for Protein ext	raction]	100ml [R0246
This product is sup cells. Proteins can b be used directly for protease inhibitors.	plied as a ready-to-use s be extracted by adding thi further analysis such as Please add a protease inh	olution for the lysis o s buffer <mark>[R0246]</mark> to th western blotting. This ibitor cocktail, if nece	of the cultured mammaliar ne cells and the extract car s product does not include essary.
Application			
Add the following pro Leupeptin 10 Pepstatin A 1 Aprotinin 3 AFBSF 1	tease inhibitors to RIPA buffer) µg/mL µg/mL } µg/mL mM	[R0246].	Add Further experiment
 Wash the cultured twice with PBS. Remove PBS and [R0246] containing manufacturer's RII inhibitors to 1.0 x Incubate the cells Centrifuge the cell Measure the prote Analyze the super 	mouse myeloma-derived cell s add 200 μ L of either cold RIPA g protease inhibitors or the othe PA buffer containing the same p 10 ⁶ cells. for 15 minutes on ice. s at 10,000 x g for 10 minutes in concentration of the superna- natants using western blotting.	sp2/0 h buffer er protease at 4 °C. atants. Washed cells	
Extracted Pr	otein Concentration	Western Blotting	
800		The extracts were tr membrane after ele Anti-β actin antibody Equal or better dete	ransferred to a PVDF ctrophoresis. y was used for detection. ction was observed than
400		that of the other ma	nutacturer's product.
400		R0246	other manufacturer's product

Y0021 is a ready-to-use solution for protein extraction from cultured *Escherichia coli* (*E. coli*) / yeast cells. By suspending cells in Y0021 and then centrifuging, the supernatant containing proteins can be obtained. Extracted protein can be used in downstream applications such as electrophoresis and western blotting.

Nervous Tissue Protein Extraction Buffer

B6279 is a ready-to-use solution for protein extraction from nervous tissue. By suspending tissue in **B6279** and then centrifuging, the supernatant containing proteins can be obtained. Extracted protein can be used in downstream applications such as electrophoresis and western blotting

100mL [B6279]

Peroxidase Substrates

TMB [for ELISA] (Ready-to-use solution) (= 3,3',5,5'-Tetramethylbenzidine (Ready-to-use solution))

100mL [T3854]

Application

- 1) Add 100µL of TMB solution [T3854] to each well.
- 2) Incubate the plate at room temperature for 30 minutes.
- 3) Add 100µL of 1N HCl solution [H1202] to each well to terminate the reaction.
- 4) Measure the absorbance of each well at 450 nm.

When **T3854** reacts with horseradish peroxidase (HRP), a blue colored soluble reaction product appears thus it can be used for ELISA.

This product cannot be used for Western blotting which needs a precipitate.



Figure. An example of use by the above method

TMB [for Western blotting] (Ready-to-use solution) (= 3,3',5,5'-Tetramethylbenzidine (Ready-to-use solution))

100mL [T3855]

Application

- 1) Incubate a blotting membrane with an HRP-conjugated antibody and then wash the membrane.
- 2) Incubate the washed membrane with TMB solution [T3855] until color development.
- 3) Add deionized water to stop color development.

When **T3855** reacts with HRP, a blue-purple precipitate appears thus it can be used for Western blotting.

This product cannot be used for ELISA which needs a soluble reaction product.



Figure. An example of Western blotting by the above method

M : molecular weight marker 1 : Target protein A 4-Chloro-1-naphthol (Ready-to-use solution) [for Western blotting] 100mL [C3384] (= 4-CN (Ready-to-use solution)) Application 1) Incubate a blotting membrane with an HRP-conjugated antibody and then wash the membrane. 2) Incubate the washed membrane with 4-CN solution [C3384] until color development. 3) Add deionized water to stop color development. kDa 70 Figure. An example of Western blotting by the above method 45 M : molecular weight marker 1 : Target protein B (Middle concentration) 2 : Target protein B (Low concentration) 32 Μ 2 1 AzBTS (Ready-to-use solution) [for ELISA] (= 2,2'-Azinobis(3-ethylbenzothiazoline-6-sulfonic Acid Ammonium Salt) (Ready-to-use solution)) 100mL [A3176]

Application

- 1) Add 100µL of AzBTS solution [A3176] to each well.
- 2) Incubate the plate at room temperature for 30 minutes.
- 3) Within 1 hour from start the reaction, measure the absorbance of each well at 405 nm.



Figure. An example of use by the above method

Related Products

Sodium Hydroxide (1mol/L in Water)	500mL [S0542]
Hydrochloric Acid (1mol/L)	500mL [H1202]
Peroxidase from Horseradish	100mg / 1g [P0073]
Horseradish Peroxidase Maleimide Conjugate (0.5mg×3)	1set [H1621]
Horseradish Peroxidase NHS Ester Conjugate (0.2mg×3)	1set [H1746]
Anti-6xHis Monoclonal Antibody (6A12) HRP Conjugate	0.05mg/1vial [A3075]
Anti-Protein A Chicken Polyclonal Antibody HRP Conjugate	0.05mg/1vial [A3187]
Anti-αGal Chicken Polyclonal Antibody HRP Conjugate	0.05mg/1vial [A3195]
Anti-NeuGc Polyclonal Antibody HRP Conjugate	0.05mg/1vial [A3397]
Goat Anti-Mouse IgG HRP Conjugate	0.1mg/1vial [G0407]
Goat Anti-Mouse IgM HRP Conjugate	0.1mg/1vial [G0417]
Goat Anti-Rabbit IgG HRP Conjugate	0.1mg/1vial [G0418]
Sheep Anti-Chicken IgY HRP Conjugate	0.1mg/1vial [S0999]
Protein A HRP Conjugate	0.2mg/1vial [P2466]
Streptavidin HRP Conjugate	0.1mg/1vial [S0972]

Soluble Substrates (for ELISA etc.)

For such as ELISA, substrates generating soluble dyes with peroxidase.

AzBTS (= 2,2'-Azinobis(3-ethylbenzothiazoline-6-sulfonic Acid Ammonium Salt))	1g <mark>[A2166]</mark>
OPD·2HCI (= 1,2-Phenylenediamine Dihydrochloride)	1g <mark>[P1144]</mark>
OPD (= 1,2-Phenylenediamine)	1g / 5g <mark>[P1805]</mark>
TMB (= 3,3',5,5'-Tetramethylbenzidine)	1g / 5g [T2573]

Soluble Substrates (for determining H_2O_2)

Substrates generating soluble dyes for determining hydrogen peroxidase (H_2O_2) by various enzyme reactions.

4-AA•2HCI (= 4-Aminoantipyrine Hydrochloride)	5g / 25g <mark>[A0257]</mark>
4-AA (= 4-Aminoantipyrine)	1g / 5g <mark>[A2254]</mark>
5-ASA (= 5-Aminosalicylic Acid) ^{*1}	5g / 25g <mark>[A2291]</mark>
DCHBS (= 3,5-Dichloro-2-hydroxybenzenesulfonic Acid Sodium Salt) *1	25g <mark>[D1928]</mark>
2,4-DCP (= 2,4-Dichlorophenol) ^{*1}	1g / 5g <mark>[D3865]</mark>
DMA (= <i>N</i> , <i>N</i> -Dimethylaniline) ^{*1}	1g / 5g <mark>[D3866]</mark>
DMT (= <i>N</i> , <i>N</i> -Diethyl- <i>m</i> -toluidine) ^{*1}	1g / 5g <mark>[D3868]</mark>
TOOS Hydrate (= Sodium 3-[Ethyl(<i>m</i> -tolyl)amino]-2-hydroxy-1-propanesulfonate Hy	'drate) ^{*1}
	1g / 5g <mark>[S0805</mark>]
ALPS (= Sodium 3-(<i>N</i> -Ethylanilino)propanesulfonate) *1	200mg / 1g <mark>[S0817]</mark>
ADOS (= Sodium 3-(<i>N</i> -Ethyl-3-methoxyanilino)-2-hydroxy-1-propanesulfonate) ^{*1}	200mg / 1g <mark>[S0826]</mark>
HDAOS (= N-(2-Hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline Sodium Salte) *1	200mg / 1g <mark>[\$0827]</mark>
MBTH•HCI (= 3-Methyl-2-benzothiazolinonehydrazone Hydrochloride)	1g / 5g [M2155]

*1 : Used together with A2254 (or A0257)

Precipitate Substrates

For such as immunohistochemical staining or immunoblotting, substrates arising precipitate products with peroxidase.

AEC (= 3-Amino-9-ethylcarbazole)	1g / 5g [A2167]
4-CN (= 4-Chloro-1-naphthol)	1g / 5g <mark>[C2291]</mark>
DAB (= 3,3'-Diaminobenzidine)	1g / 5g <mark>[D3756]</mark>
DAB•4HCI (= 3,3'-Diaminobenzidine Tetrahydrochloride Hydrate)	1g / 5g <mark>[D3757]</mark>
o-Dianisidine *2	1g / 5g <mark>[D3864]</mark>
o-Dianisidine Dihydrochloride *2	1g / 5g <mark>[D3893]</mark>
DMPD·2HCI (= <i>N</i> , <i>N</i> -Dimethyl-1,4-phenylenediamine Dihydrochloride) *3	1g / 5g <mark>[D3931]</mark>
1-Naphthol * ³	1g / 5g [N0864]
*2 : By combinating N0864 and D3931 *3 : Used together with C2291	



Soluble Substrates

4-Nitrophenyl Phosphate Disodium Salt Hexahydrate	1g / 5g [D4005]
4-Nitrophenyl Phosphate Di(tris) Salt Hydrate	5g / 25g [N0422]
1-Naphthylphosphoric Acid Monosodium Salt Monohydrate	1g / 5g / 25g [N0452]
1-Naphthylphosphoric Acid Disodium Salt Hydrate	1g / 5g [P0263]

Precipitate Substrates

For such as immunohistochemical staining or immunoblotting, substrates arising precipitate dyes with alkaline phosphatase.

5g / 25g	[B0785]
100mg / 1g	[B1239]
1g / 5g	[B3581]
200mg	[C2250]
100mg / 1g	[D0844]
100mg / 1g	[10781]
100mg / 1g	[T0250]
	5g / 25g 100mg / 1g 1g / 5g 200mg 100mg / 1g 100mg / 1g 100mg / 1g

Protein Staining Reagent

Coomassie Brilliant Blue G-250 (Ready-to-use solution) [for Electrophoresis]

500mL [C3488]

Application

- 1) After electrophoresis, wash the gel with deionized water for 5 minutes three times.
- 2) Remove the water used for washing, add C3488 till the gel is soaked, and stain the gel for 1 hour while shaking gently at room temperature.
- 3) Remove the staining solution, destain the gel with deionized water for 1 hour and check it.
- 4) If the background is high, destain the gel with deionized water overnight at room temperature.



Figure. Proteins stained by the above method (destained overnight)

Related Products

2X SDS-PAGE Sample Buffer (2-Mercaptoethanol free)	25mL [B5834]
4X SDS-PAGE Sample Buffer (2-Mercaptoethanol free)	20mL [B6104]
6X Sample Buffer (2-Mercaptoethanol free)	10mL [B6105]
Pyrogallol Red [for Protein Research]	1g [P1976]
Streptomycin Sulfate [for Protein Research]	5g/25g <mark>[\$0834]</mark>
Acrylamide Monomer [for Electrophoresis]	25g / 500g [A1132]
30% Acrylamide / Bis-acrylamide (29:1) [for Electrophoresis]	250mL [A3217]
30% Acrylamide / Bis-acrylamide (37.5:1) [for Electrophoresis]	250mL [A3218]
Acid Black 1 [for Electrophoresis]	5g [A2097]
Ammonium Peroxodisulfate [for Protein Research]	5g / 25g <mark>[A2098]</mark>
Coomassie Brilliant Blue G-250 [for Electrophoresis]	5g [B3193]
Coomassie Brilliant Blue R-250 [for Electrophoresis]	5g [B3194]
Bromophenol Blue Sodium Salt [for Electrophoresis]	1g [B3195]
Sodium Deoxycholate [for Electrophoresis]	25g [D1820]
DL-Dithiothreitol [for Electrophoresis]	1g/5g [D3647]
Glycerol [for Electrophoresis]	1g [G0316]
Glycine [for Electrophoresis]	25g / 500g [G0317]
N,N'-Methylenebisacrylamide [for Electrophoresis]	25g / 100g [M0506]
2-Mercaptoethanol [for Electrophoresis]	5g / 25g [<mark>M1948]</mark>
Sodium Dodecyl Sulfate (=SDS) [for Electrophoresis]	25g / 500g [S0588]
N,N,N',N'-Tetramethylethylenediamine (=TEMED) [for Electrophoresis]	5g / 25g [T2515]
Tris(hydroxymethyl)aminomethane (=Tris-Base) [for Electrophoresis]	25g/500g [T2516]

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волгоград (844)278-03-48 Вологда (8172)26-4159 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89 Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48

Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81 Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73

Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47 Ростов-на-Дону (863)308-18-15 Рязань (4912)46-61-64 Самара (846)206-03-16 Саранск (8342)22-96-24 Санкт-Петербург (812)309-46-40 Саратов (845)249-38-78 Севастополь (8692)22-31-93 Симферополь (3652)67-13-56 Смоленск (4812)29-41-54 Сочи (862)225-72-31 Ставрополь (8652)20-65-13 Сургут (3462)77-98-35

Сыктывкар (8212)25-95-17 Тамбов (4752)50-40-97 Тверь (4822)63-31-35 Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удэ (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

Адрес: https://tci.nt-rt.ru/ || эл.почта: tic@nt-rt.ru