

Реагент для очистки тканей

ЖИВОТНЫХ

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-4159
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48

Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81
Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73

Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35

Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

Адрес: <https://tci.nt-rt.ru/> || эл.почта: tic@nt-rt.ru

Animal Tissue-Clearing Reagent CUBIC

Products

CUBIC trial kit (including mounting solution)

; All-in-One basic tissue clearing kit containing CUBIC-L (50mL), CUBIC-R+(M) (50mL), and RI-Matched Mounting Solution (RI = 1.520, 50mL) 1kit [C3942]

CUBIC-L (for delipidation and decoloring) 25mL / 100mL / 500mL [T3740]

CUBIC-R+(N) (for RI matching) 25mL / 100mL / 500mL [T3983]

CUBIC-R+(M) (for RI matching) 25mL / 100mL [T3741]

CUBIC-B (for decalcification) 25mL / 100mL [T3780]

CUBIC-HL (for delipidation strongly and quenching autofluorescence) 25mL / 100mL [T3781]

CUBIC-P (with perfusion before tissue excision) 25mL / 100mL [T3782]

CUBIC-X1 (for expansion) 25mL / 100mL [T3866]

CUBIC-X2 (for RI matching with expansion) 25mL / 100mL [T3867]

CUBIC-HV™1 3D immunostaining kit (Casein separately) 1kit [C3717]

CUBIC-HV™1 3D nuclear staining kit 1kit [C3709]

Related Products

Mounting Solution (RI 1.520) [for CUBIC-R+] 50mL [M3294]


Mounting Solution (RI 1.467) [for CUBIC-X2] 50mL [M3292]

Mouse Anti-NeuN Monoclonal Antibody 0.1mg/vial [M3586]

Goat Anti-Mouse IgG₁ Fab Fragment Cyanine 3 Conjugate 0.05mg/vial [G0598]



Whole-brain clearing



Whole-body clearing
with nuclei staining and immunostaining

These products were developed by Prof. Hiroki R. Ueda (The University of Tokyo / RIKEN) and are under invention licenses by RIKEN, Japan.

*CUBIC-HV™ is a registered trademark of CUBICStars Co.

Advantages

- **Basic protocol ;**

Clearing of whole mouse bodies as well as animal organs can be achieved by using two reagents in sequence: CUBIC-L [T3740] for delipidation and either CUBIC-R+(N) [T3983] or CUBIC-R+(M) [T3741] for RI matching.

The difference between CUBIC-R+(N) [T3983] and CUBIC-R+(M) [T3741]:

CUBIC-R+(N) is inexpensive and easier to handle because it raises less precipitation.

The fluorescence signal may decay, but the fluorescence signals of samples in CUBIC-R+(N) can be observed for several days after immersion. CUBIC-R+(M) is superior in retaining the fluorescence signal. However, at low temperatures such as in winter, it may precipitate. In that case, it can be resolved by placing the sample at 37°C for a few days.

For these reasons, it is recommended to try CUBIC-R+(N) first and then use CUBIC-R+(M) if fluorescence signal cannot be found.

- **Optional protocol ;**

The following products can easily clear tissues, such as bones or highly fatty tissues which were previously difficult to clear.

CUBIC-B [T3780] for bone

CUBIC-HL [T3781] for highly fatty tissues

- **For efficiently aiding with perfusion fixation for mouse perfusion ;**

CUBIC-P [T3782]

- **Expansion protocol ;**

The following products can clear tissues with expansion.

CUBIC-X1 [T3866] for expansion tissues

CUBIC-X2 [T3867] for RI matching with keeping the expanded size of tissues

- **For staining thick and large specimens uniformly ;**


CUBIC-HV™ 1 3D immunostaining kit [C3717] for 3D immunostaining

CUBIC-HV™ 1 3D nuclear staining kit [C3709] for 3D nuclear staining

- **Tissue expansion enables acquisition of images easy.**

- **Preserve the fluorescent protein signals except CUBIC-HL [T3781].**

- **Using light-sheet fluorescent microscopy (LSFM) or confocal laser-scanning microscopy (CLSM) enables the whole-organ / body imaging at a cellular resolution.**



Whole-body clearing



Whole-body clearing
with propidium iodide staining

Direction for Use : Mouse whole-organ clearing protocol

Fix 4% PFA 1 day	Wash x 3 PBS > 2 hr x 3	Delipidation 50% CUBIC-L 6 – 24 hr	Delipidation CUBIC-L > 2 days	Wash x 3 PBS > 2 hr x 3	(Staining) Stains > 3 days	(Wash x 3) PBS > 2 hr x 3	(Fixation) 1% FA 1 day	(Fixation) 1% FA 1 hr	(Wash x 3) PBS > 2 hr x 3	RI match 50% CUBIC-R+ 1 day	RI match CUBIC-R+ > 1 day
------------------------	-------------------------------	--	-------------------------------------	-------------------------------	----------------------------------	---------------------------------	------------------------------	-----------------------------	---------------------------------	-----------------------------------	---------------------------------

Process	Reagent	Temp.	Time	Notes
Tissue excision				After perfusion fixation
Tissue Fix	4% PFA in PBS	4°C	1 day	
Wash x 3	PBS	RT	> 2 hr x 3	Shake gently (Same in following steps). Total 1 day
(Delipidation)	50% CUBIC-L	37°C or RT	6 – 24 hr	1:1 mixture of water and CUBIC-L. Optional
Delipidation	CUBIC-L	37°C	> 2 days	Refresh CUBIC-L on day 1, day 2 and every 2 days after day 4.
Wash x 3	PBS	RT	> 2 hr x 3	Total 1 day
(Staining)	Stains	RT	> 3 days	Optional
(Wash x 3)	PBS	RT	> 2 hr x 3	Total 1 day, When stained
(Fixation)	1% FA	4°C	1 day	Diluted of 37% FA by PBS. When stained.
(Fixation)	1% FA	37°C	1 hr	When stained.
(Wash x 3)	PBS	RT	> 2 hr x 3	When stained.
RI match	*50% CUBIC-R+	RT	1 day	1:1 mixture of water and CUBIC-R+
RI match	*CUBIC-R+	RT	> 1 day	

*Both CUBIC-R+(N) [T3983] and CUBIC-R+(M) [T3741] can be used.

Application

● An adult mouse brain after excision



● After pre-treatment step of 4 mL 50% CUBIC-L at room temperature overnight



● After delipidation step of 4 mL CUBIC-L at 37°C for 5 days (Refresh CUBIC-L on day 1, day 2 and day 4)



● After pre-treatment step of 4 mL CUBIC-R+(M) at room temperature overnight



● Observation the sample in Mounting Solution (RI = 1.520) [M3294] after RI matching of 4 mL CUBIC-R+(M) at room temperature overnight



Each sample of these images is immersed in each reagent.

- Light penetrates the organ.
 - CUBIC-L does not get colored after treatment.
- Above points are the signs of end of delipidation.

Total

- CUBIC-L : 14 mL
- CUBIC-R+(M) : 6 mL

The reagent volumes of the left example is in the case of usage in a 5 mL-tube.

Work in a tube whose diameter is a little larger than that of organs and the volume of reagents is half of that of tubes.

PFA : paraformaldehyde, RT : room temperature

Animal Tissue-Clearing Reagent CUBIC

Direction for Use : Mouse whole-brain clearing with expansion protocol

Fix 4% PFA 1 day	Wash x 3 PBS > 2 hr x 3	Delipidation 50% CUBIC-L 3 hr	Delipidation CUBIC-L 5 - 14 days	Wash PBS 1 day	Staining Stains 3 days	Wash PBS 1 day	Fixation 1% FA 1 day	Fixation 1% FA 1 hr	Wash x 3 PBS > 2 hr x 3	Swelling CUBIC-X1 2.5 days	RI match CUBIC-X2 1.5 days
------------------------	-------------------------------	-------------------------------------	--	----------------------	------------------------------	----------------------	----------------------------	---------------------------	-------------------------------	----------------------------------	----------------------------------

Process	Reagent	Temp.	Time	Notes
Tissue excision	After perfusion fixation			
Tissue Fix	4% PFA in PBS	4°C	1 day	
Wash x 3	PBS	RT	> 2 hr x 3	Shake gently (Same in following steps). Total 1 day
Delipidation	50% CUBIC-L	37°C	3 hr	1:1 mixture of water and CUBIC-L.
Delipidation	CUBIC-L	37°C	5 - 14 days	Refresh CUBIC-L every 4 days. 5 days for 1-week-old mice 7 days for 3-week-old mice 14 days for 8-week-old and 6-month-old mice
Wash	PBS	RT	1 day	
Staining	Stains	RT	3 days	
Wash	PBS	RT	1 day	
Fixation	1% FA	4°C	1 day	Diluted of 37% FA by PBS.
Fixation	1% FA	37°C	1 hr	
Wash x 3	PBS	RT	> 2 hr x 3	
Swelling	CUBIC-X1	4°C	2.5 days	
RI match	CUBIC-X2	RT	1.5 days	Refresh CUBIC-X2 every 12 hours.

Application

- Example of usage of mouse brain clearing and expansion.
- Pre-treatment step of 3 mL 50% CUBIC-L at 37°C for 3 hours after PBS wash.
- Delipidation step of 3 mL CUBIC-L at 37°C for 14 days.
(Refresh CUBIC-L on day 4, day 8 and day 12.)
- Wash by PBS, staining by staining reagents and wash by PBS.
- Expansion step of 30 mL CUBIC-X1 at 4°C for 2.5 days.
- Observation the sample in Mounting Solution (RI = 1.467) [M3292] after RI matching of 40 mL CUBIC-X2 at room temperature for 1.5 days.
(Refresh CUBIC-X2 every 12 hours.)


Total

- CUBIC-L : 10.5 mL
- CUBIC-X1 : 30 mL
- CUBIC-X2 : 120 mL

*For nuclear staining, use 30 µg/mL Propidium iodide (PI) and 1.5 M NaCl in PBS.

Since the expanded brains are fragile, careful handling is required after the swelling step.

PFA : paraformaldehyde, RT : room temperature



Magnified view of a transgenic mouse brain after clearing-expansion protocol

3D Tissue Staining

3D Tissue Staining Kits CUBIC-HV™

Introduction

- Stain bulky specimens uniformly.
(Includes two nuclear stains and an antibody control)
- CUBIC-L [T3740] and either CUBIC-R+(N) [T3983] or CUBIC-R+(M) [T3741] (sold separately) required for upstream / downstream sample processing.



Kit Components

CUBIC-HV™1 3D immunostaining kit (Casein separately)

1kit [C3717]

- 2 x Immunostaining Buffer (Casein separately) (for 10 tests)
- 1 x Immunostaining Washing Buffer (for 10 tests)
- 10 x Immunostaining Additive (for 10 tests)
- Anti NeuN Mouse IgG1 Antibody (1mg/mL) (for 2 tests)
- Subdivided Casein (1 vial)
- 10 packs of 15mL tube

CUBIC-HV™1 3D nuclear staining kit

1kit [C3709]

- 1 x 3D Nuclear Staining Buffer (for 10 tests)
- 100 x 3D Nuclear Staining Washing Buffer (for 10 tests)
- 200 x DAPI 2HCl (1mg/mL in Water) [for Cell Staining] (for 10 tests)
- 100 x Propidium Iodide (1mg/mL in Water) [for Cell Staining] (for 10 tests)
- 10 packs of 5mL tube

These volumes are for mouse adult brains. Components of the kits are subject to change without notice.

По вопросам продаж и поддержки обращайтесь:

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-4159
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48

Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81
Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73

Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35

Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

Адрес: <https://tci.nt-rt.ru/> || эл.почта: tic@nt-rt.ru